

BagMixer® Lab blenders range

BagMixer® are highly efficient lab blenders for solid sample preparation. With a sterile bag, they guarantee optimal microbiological extraction in less than 60 seconds and eliminate all risks of cross contamination.

More than **15,000** BagMixer® used daily!



MiniMix® 100
5-80 mL
Sample up to 10 g




BagMixer® 400
50-400 mL
Sample up to 40 g



JumboMix® 3500
200-3500 mL
Sample up to 400 g

- ▶ World best selling lab blender
- ▶ 11 models for most applications (from 5 mL to 3500 mL)
- ▶ 35 years experience
- ▶ Life-time service and spare parts



MiniMix® 100 100 mL lab blenders

The MiniMix® 100 is designed for the blending of small samples. Thanks to the adjustable paddles, your blending is more accurate and efficient.



MiniMix® 100 P CC

- Multifunction digital screen
- Variable speed (3, 6, 9, 12 strokes/sec.)
- Variable time (1 sec. - 1 h or ∞)
- Adjustable blending power (patented)
- Q-Tight® closing
- Click & Clean® removable paddles (patented)
- Side-by-side paddle stop

MiniMix® 100 W CC

- Multifunction digital screen
- Variable speed (3, 6, 9, 12 strokes/sec.)
- Variable time (1 sec. - 1 h or ∞)
- Adjustable blending power (patented)
- Q-Tight® closing
- Click & Clean® removable paddles (patented)
- Side-by-side paddle stop
- **Window door**
- **Security DripTray**

Discover the MiniMix® accessories p. 22!

Delivered with: Power cord, user's manual, security DripTray®, free sterile packs of filter bags. (*only on MiniMix® 100 W CC)

Technical specifications

MiniMix® 100	P CC	W CC
Useful volume	5 - 80 mL	
Power	110/220V~ 50-60Hz	
Size (w x d x h)	20 x 37.5 x 20 cm	
Weight	12 kg	11.8 kg
Reference	011 230	013 230

Ideal with the sterile bags



BagPage®
80 mL



BagPage®
100 mL



BagLight®
100 mL

Applications

- Homogenization of all types of samples up to 10 g
- Clinical, pharmaceutical
- Veterinary (biopsy, organ tissue sampling...)
- Molecular extraction (proteins, lipids, DNA...)
- Toxin research (HPLC, GCMS...)

Advantages


- ▶ Silent and resistant
- ▶ Adjustable blending power
- ▶ Easy cleaning: door opens at 270°

BagMixer[®] 400 Unique features by interscience

Adjustable blending power

- Up to 28 kg of pressure per paddle

PATENTED



Small/hard samples:
Move the paddles forward
for more pressure



Big/soft samples:
Move the paddles backwards
for less pressure

Super Silent Blending⁽¹⁾

- Less than 48 dB
- Double layer window door



ISO 3744 acoustic test

Q-Tight[®] system

- Applies 250 kg of pressure to avoid leaks from the top

Easy-to-clean

- TotalAccess System: open the door flat on the bench⁽¹⁾
 - Click & Clean[®] removable paddles
- PATENTED**

Robust and reliable

- High quality 304L stainless steel
- 3 year guarantee
- Lifetime guarantee on window and shock absorbers

Advanced lighting system⁽¹⁾


- LightCode System: indicates the state of the blending process
- LED lighting inside the blending chamber

Ultra stable⁽¹⁾

- Low vibration
- Can be used next to a dilutor

LiquidSensor⁽²⁾

- Automatic stop of the blending in case of bag leakage
- PATENTED**



(1) Only on BagMixer[®] S&SW (2) Only on BagMixer[®] SW

BagMixer® 400 400 mL lab blenders

World best-seller in labs, BagMixer® 400 is an efficient lab blender for all microbiological analyses.

ISO 7218	ISO 6887-1	FDA BAM Bacteriological Analytical Manual
3 YEAR GUARANTEE	SHOCK ABSORBERS LIFETIME GUARANTEE	WINDOW LIFETIME GUARANTEE



Less than
48 dB!



Less than
48 dB!

BagMixer® 400 P

- Fixed speed (8 strokes/sec.)
- Variable time (30-210 secs. or ∞)
- Q-Tight® closing
- Also available in CoolMix® version

BagMixer® 400 W

- Fixed speed (8 strokes/sec.)
- Variable time (30-210 secs. or ∞)
- Q-Tight® closing
- **Window door**
- **Security DripTray**

BagMixer® 400 CC

- Window door
- Security DripTray
- Q-Tight® closing
- **Multifunction digital screen**
- **Variable speed (4, 6, 8, 10 strokes/sec.)**
- **Variable time (1 s - 1 h or ∞/count-down)**
- **Adjustable blending power (patented)**
- **Side-by-side paddle stop**
- **Silent brushless motor**
- **Click & Clean® removable paddles (patented)**

BagMixer® 400 S

- Multifunction digital screen
- Variable speed (4, 6, 8, 10 strokes/sec.)
- Variable time (1 s - 1 h or ∞/count-down)
- Adjustable blending power (patented)
- Q-Tight® closing
- Side-by-side paddle stop
- Click & Clean® removable paddles (patented)
- **Super silent brushless motor**
- **Integrated security DripTray**
- **TotalAccess System**
- **LightCode System**

BagMixer® 400 SW

- Multifunction digital screen
- Variable speed (4, 6, 8, 10 strokes/sec.)
- Variable time (1 s - 1 h or ∞/count-down)
- Adjustable blending power (patented)
- Q-Tight® closing
- Side-by-side paddle stop
- Click & Clean® removable paddles (patented)
- Super silent brushless motor
- Integrated security DripTray
- TotalAccess System
- LightCode System
- **Double layer window door**
- **Liquid Sensor**

Demo video
www.interscience.com

Discover the BagMixer® accessories p. 22!

Delivered with: Power cord, user's manual, security DripTray®, free sterile packs of filter bags. (*in option for BagMixer® 400 P)

Technical specifications

BagMixer® 400	P	W	CC	S	SW
Useful volume	50 - 400 mL				
Power	110-127/220-240 V~ 50-60Hz		100-240 V~ 50-60Hz		
Size (w x d x h)	27 x 38.5 x 28 cm		27 x 41.5 x 26.5 cm	27.5 x 46.5 x 36 cm	
Weight	16.6 kg		16 kg	25.2 kg	25 kg
Reference	021 230	022 230	024 230	025 000	025 100

Ideal with the sterile bags



Applications

- Blending of all types of samples up to 40 g
- Microbiology: food, pharmaceutical, clinical, veterinary, environmental...

Advantages

- ▶ **Silent and robust**
- ▶ **Adjustable blending power**
- ▶ **Easy cleaning**

JumboMix® 3500 3500 mL lab blenders

The JumboMix® 3500 enables fast blending of large samples for numerous applications in microbiology.



JumboMix® 3500 VP

- Multifunction digital screen
- Variable speed (1.5, 3, 4.5, 6 strokes/sec.)
- Variable time (1 sec. - 1 h or ∞)
- Adjustable blending power (patented)
- Q-Tight® closing
- Side-by-side paddle stop



JumboMix® 3500 VW

- Multifunction digital screen
- Variable speed (1.5, 3, 4.5, 6 strokes/sec.)
- Variable time (1 sec. - 1 h or ∞)
- Adjustable blending power (patented)
- Q-Tight® closing
- Side-by-side paddle stop
- **Window door**
- **Security DripTray**



Special Trichinella

JumboMix® 3500 WarmMix®

- Multifunction digital screen
- Variable speed (1.5, 3, 4.5, 6 strokes/sec.)
- Variable time (1 sec. - 1 h or ∞)
- Adjustable blending power (patented)
- Q-Tight® closing
- Side-by-side paddle stop
- Security DripTray
- **Heated door (20°C to 50°C)**

Demo video www.interscience.com

Applications

- 375 g for salmonella
- 75, 325 or 375 g sample for E. Coli O157:H7
- 125 g sample for Listeria and Listeria monocytogenes
- Sample pooling (e.g.: 15 samples of 25 g to reach 375 g)

Discover the JumboMix® accessories p. 22!

Delivered with: Power cord, user's manual, security DripTray®, free sterile packs of filter bags. (*only on JumboMix® 3500 VW and JumboMix® 3500 WarmMix)

Technical specifications

JumboMix® 3500	VP	VW	WarmMix®
Useful volume	200 - 3500 mL		
Power	110/220V~ 50-60Hz		
Size (w x d x h)	51.8 x 47 x 47 cm		63 x 50 x 47 cm
Weight	53.2 kg	50.2 kg	63 kg
Reference	031 230	032 230	033 230

Ideal with the sterile bags



Option



Advantages

- ▶ Silent & robust
- ▶ Adjustable blending power
- ▶ Easy cleaning: door opens at 270°